

Can a Mixing Inlet Interposed Between a Pressurized Metered Dose Inhaler-Spacer Combination and Cascade Impactor Enable Effective Performance Testing at Widely Different Inhalation Flow Rates?

Mark W. Nagel¹, Jason A. Suggett¹, Rubina Ali¹, and Jolyon P. Mitchell²
¹ Trudell Medical International, London, Ontario, Canada; ² Jolyon Mitchell Inhaler Consulting Services Inc., London, Ontario, Canada

SUMMARY

The MI-NGI combination was in good agreement with the NGI at 30 L/min for almost all measures.

CLINICAL BACKGROUND

- Multi-stage cascade impactors (CIs) are used to characterize of emitted aerosol aerodynamic particle size distribution (APSD) from orally inhaled products (OIPs).^{1,2}
- The Next Generation Impactor (NGI) [3] operates at 15 L/min [4] and between 30 and 100 L/min⁵
- The low flow Marple-Miller impactor (MMI) was developed primarily to enable low flow rate characterizations at 4.9 and 12.0 L/min in support of delivery devices for infants and small children.⁶
- The Miller design of mixing inlet (MI)⁷ is located between the mouthpiece of the OIP-on-test connected to the induction port and the NGI.
 - The MI allows the CI to operate at an optimum flow rate for the inhaler on test.
 - The flow rate leaving the inhaler can be reduced for pediatric device testing by adding a fixed flow rate of clean make-up air less than that required by the CI via its sidearm.⁸

STUDY PURPOSE

The primary purpose of our study was to discover whether a Miller MI used in conjunction with an NGI could be used to evaluate a pMDI with spacer add-on device at 30 L/min and at sampling flow rates comparable with those for the low flow MMI.

MATERIALS AND METHODS

- A purpose-built tube spacer was used in conjunction with the albuterol pMDI as a common aerosol source to enable equal comparison of the different equipment configurations to be made with an add-on device present.
- Ventolin* (GSK, Canada) pMDIs delivering 100 µg albuterol base equivalent ex valve per actuation were fitted to the spacer exit (ca. 149 ml) as shown in Figure 1.
- Replicate measurements (n = 5) were made at each sampling condition
- We first made baseline APSD measurements with the spacer mouthpiece attached directly to the NGI via a USP/PhEur induction port, sampling from the inhaler-spacer at 30 L/min determined by mass flowmeter (model 4040, TSI Inc., St. Paul, MN).
- We then as sampled at 4.9 and 12 L/min via the low flow MMI connected directly to the spacer mouthpiece via a USP/PhEur induction port
- Finally, with the MI interposed between the inhaler-spacer and NGI whose flow rate was fixed at 60 L/min, we sampled from the spacer at 4.9, 12, and 30 L/min in separate measurements.
- A supply of compressed air was fed to the MI sidearm so that the nominal flow rate measured by a mass flowmeter located at the entry to the MI could be adjusted to the desired sampling flow rate
- The mass of albuterol recovered from the induction port and each Brij-35 surfactant-coated cup of the NGI or low flow MMI was subsequently determined by a validated UV-HPLC method.




Figure 1: pMDI-Tube Spacer Configuration Shown with MI-NGI Configuration

DISCUSSION

- Mass balances at all sampling conditions were close to the expected value of 100 µg albuterol/actuation, indicating minimal losses had occurred within each of the sampling arrangements.
 - TM_{exHC}** sampling via the low flow MMI at 4.9 L/min (27.2 ± 5.1 µg/actuation) and 12 L/min (40.5 ± 8.2 µg/actuation) were lower than equivalent values using the NGI-MI combination (32.3 ± 4.8 and 50.9 ± 5.0 µg/actuation) [2-tailed, un-paired t-test at each flow rate, p ≥ 0.04].
 - Larger internal losses in the low flow MMI may explain the divergence.
 - TM_{exHC}** (MI-NGI combination, 53.7 ± 3.5 µg/actuation) agreed with the corresponding values using the NGI alone (56.0 ± 3.9 µg/actuation) [un-paired t-test, p = 0.35], confirming that losses in the MI were minimal.
- FPM_{<4.7µm} sampling via the low flow MMI at either 4.9 L/min (22.3 ± 4.3 µg/actuation) or 12 L/min (34.4 ± 7.8 µg/actuation), were equivalent to FPM_{<4.5µm} using the MI-NGI (27.7 ± 3.8 (4.9 L/min) and 43.4 ± 4.3 (12 L/min) µg/actuation respectively) [p ≥ 0.054].
- Low-flow MMI-derived values of FPF_{<4.7µm} at 4.9 L/min (91.9 ± 1.7%) were close to FPF_{<4.5µm} using the MI-NGI combination at 12 L/min (91.2 ± 1.4%) [p = 0.31].
 - The cause of the small divergence observed at 12 L/min (88.9 ± 0.8 low-flow MMI), compared with the NGI-MI combination (92.4 ± 0.6%) [p < 0.001], may be linked to the lower upper limit for fine particles in the low flow MMI data.
- FPM_{<4.5 µm} at 30 L/min using the MI-NGI (43.9 ± 3.3 µg/actuation) was comparable with FPM_{<4.0 µm} obtained by the NGI alone (45.9 ± 4.1 µg/actuation), consistent with the findings for **TM_{exHC}** (p = 0.42).
 - Likewise, FPF_{<4.5 µm} (87.6 ± 0.7%) determined by MI-NGI concurred with FPF_{<4.0 µm} obtained directly by the NGI (85.9 ± 2.6%) [p = 0.14].
- Importantly, estimated values of MMAD by interpolation of FPF values for the two stages either side of the median FPF associated with all the MI-NGI measurements ranged from 2.5 to 2.7 µm.
 - This range was consistent with MMAD values from the low-flow MMI.
- Spread factors determined for all sampling conditions were within 1.7 ± 0.05.

CONCLUSIONS

- The MI-NGI combination was in good agreement with the NGI at 30 L/min for almost all measures.
- FPM values using the NGI-MI sampling either at 4.9 or 12 L/min, were slightly greater than equivalent results via the low flow MMI.
 - This outcome is likely unimportant for predicting OIP performance
- We conclude that the MI-NGI combination can be used as an alternative for the low flow NGI, as well as for measurements sampling from the inhaler at 30 L/min, with the NGI operating at 60 L/min.

RESULTS

- Total recovered mass of albuterol (n = 5 replicates/condition) are presented in Table 1.

Table 1: Mass Balances for pMDI-Spacing Tube-Delivered Albuterol (mean ± SD) Sampling via Low Flow MMI, NGI alone and NGI-MI

* Includes the mass recovered from the mouthpiece of the pMDI and the spacer interior

Sampling Condition	Low Flow MMI		MI+NGI			NGI
	4.9	12	4.9	12	30	
Flow rate through CI (L/min)	60		60			30
Entry flow rate (L/min)	60		4.9	12	30	60
Total recovered* mass (µg/actuation)	95.3 ± 4.8	102.7 ± 6.7	100.4 ± 4.7	104.4 ± 4.3	100.5 ± 1.1	101.6 ± 4.0

- CI-derived measures of total mass ex VHC (**TM_{exHC}**), fine particle mass (FPM), fine particle mass fraction (FPF) and the individual upper bound sizes for fine particles are summarized in Table 2.
- We did not interpolate FPM and FPF to a fixed upper bound size to preserve the statistical information available from the measured values at each sampling condition.
- However, the cut-point size of the CI stage for each condition was chosen to be as close as possible to 5 µm aerodynamic diameter, as specified in the European Pharmacopoeia.²
- We determined MMAD values by interpolation as the size corresponding to the 50th percentile of the cumulative mass-weighted APSDs.

Table 2: CI-Derived Measures for pMDI-Spacing Tube-Delivered Albuterol (mean ± SD) Sampling Directly via Low Flow MMI, by NGI alone and NGI-MI

Sampling Condition	Low Flow MMI		MI+NGI			NGI
	4.9	12	4.9	12	30	
Flow rate through CI (L/min)	60		60			30
Entry flow rate (L/min)	60		4.9	12	30	60
TM _{exHC} (µg/actuation)	95.3 ± 4.8	102.7 ± 6.7	100.4 ± 4.7	104.4 ± 4.3	100.5 ± 1.1	101.6 ± 4.0
Upper limit for FPM (µm)	27.2 ± 5.1	40.5 ± 8.2	32.3 ± 4.8	50.9 ± 5.0	53.7 ± 3.5	56.0 ± 3.9
FPM (µg/actuation)	4.7	4.7	4.5	4.5	4.5	4.0
FPF (%)	91.9 ± 1.7	92.4 ± 0.6	91.2 ± 1.4	88.9 ± 0.8	87.6 ± 0.7	85.9 ± 2.6
MMAD (µm)	2.6	2.7	2.5	2.5	2.5	2.4
Spread Factor	1.7	1.7	1.7	1.7	1.7	1.7

- We found that all APSDs were unimodal but deviated sufficiently from lognormality, especially for the 4.9 L/min sampling flow rate conditions.
- We therefore calculated a spread factor from the square root of the 80th to 20th mass percentiles as an alternative measure to compare APSD spread.