

On-the-go with
Aero2go[™]
Chamber

Designed with the help of
people with asthma and COPD
to suit their **active** lifestyles.



2-in-1 Chamber + Protective Case

Inhaler users involved at every stage of development



What we *learned*:


- Asthma attacks often occur outside the home when people are exposed to triggers
- Difficult to inhale medication during an exacerbation
- Debris inhaled from inhaler
- Majority of adults don't use a chamber outside of their home

What they *preferred*:

- **Chamber + Protective Case (2-in-1)** preferred to smaller collapsible chamber
- Pocket-size chamber
- Discrete
- Does not look like a medical device




Reasons for *not* using a chamber¹:



Improved adherence and high patient satisfaction²

Aero2go
Chamber

Greater than **2x increase**
in spacer use **on-the-go**²



“I like how easy it stored my puffer inside so it takes up less space in my purse... what a fantastic idea!”

Becky

“It’s more discrete and easy to carry.”

David

“Easy to clean. Easy to carry. Easy to use.”

Rose



Patient-inspired, real life-ready

Aero2go
Chamber

Easy to use

- Simply breathe normally to inhale the medication
- Inhalation valve holds medication suspended until ready to inhale

Easy to access in emergency

- Unique opening mechanism provides quick access to inhaler
- Cap secures the chamber in closed position and helps prevent it from inadvertently opening

Chamber also functions as a protective case for the inhaler


- Only need to carry 1 item versus 2
- Keeps the inhaler clean and secure
- Protects the inhaler (chamber enclosed on both ends)

Durable

- Shatter resistant
- Can be used directly out of package

Anti-static materials

Now fits even more types of metered dose inhalers




2 in 1 Chamber +
Protective Case



Proven to deliver inhaler medication better than the inhaler alone³


Aero2go
Chamber

Chambers optimize medication delivery




When the inhaler is used alone, even a one second inhalation delay reduces medication delivery by **80%**

AeroChamber2go™ + Inhaler




Provides consistent medication delivery, even after time delays



Almost **10x** more medication deposits in back of throat⁴

Inhaler alone



More medication is delivered to the lungs⁵

Inhaler with **AeroChamber®** brand of chambers

Chambers for use at home and on-the-go

Aero2go
Chamber

Part of the **AeroChamber®** Family of Chambers



REFERENCES:

1. TMI data on file. Survey of 715 inhaler patients 2018. 2. Suggett J. Improved adherence and high patient satisfaction with 'on-the-go' spacer. Canadian Respiratory Conference 2023. 3. Suggett J et al. Patient centered development of a new valved holding chamber designed specifically for on-the-go use. European Respiratory Journal 2020 56: 3181; DOI: 10.1183/13993003.congress-2020.3181 4. Suggett J *et al.* Drug Delivery to the Lungs 2017. 5. Gillen M *et al.* Pulmonary Pharmacology and Therapeutics 2018 Oct;52:7-17.



**TRUDELL MEDICAL
INTERNATIONAL**

TM/® trademarks and registered trademarks of Trudell Medical International (TMI). *Trademarks of the respective companies.
© TMI 2025. All rights reserved. Manufactured in Canada with Canadian and USA parts.
MD-518A-0125



Designed and
Developed in Canada
with the help of
people like you

Portfolio >